



Transfection reagent

M *In Vivo*
PolyMag

In Vivo Nucleic Acids delivery

Protocol



Magnetofection Technology

This reagent needs to be used with specific magnets

IMPORTANT NOTES – Before you begin

1. The conditions provided above might require some further optimizations depending on your nucleic acids, animal, territory, routes of injection etc...
2. It is suggested to use 1 μ L of *In vivo* PolyMag in initial experiments.
3. Allow reagents to reach RT and gently vortex them before forming complexes.
4. The final DNA concentration should not exceed 0.5 mg/mL.
5. Dilutes the reagents with deionized water for doses less than 1 μ L. Discard the diluted reagent after use.
6. Nucleic acids should be as pure as possible, endotoxins free and prepare in water
7. For the complexes preparation and injection, prefer glucose 5% solution or saline buffer (HBS, PBS, normal saline, Ringer's solution).
8. Do not inject more than 1 mL of *In vivo* PolyMag per animal.
9. Do not inject complexes if precipitate has formed
10. Do not freeze magnetic nanoparticles
11. Do not add anything to the stock solution of magnetic nanoparticles
12. Magnet manipulation:
 - a. Manipulate carefully the magnets. Danger of injury by strong magnetic attraction of ferromagnetic material
 - b. Keep away from electronic devices and magnetic storage devices
 - c. Persons with cardiac pacemakers should not work with these magnets

For additional information and protocols (optimization, scaling, co-transfection...) tips, troubleshooting or other applications



www.ozbiosciences.com

Any questions?



tech@ozbiosciences.com

In vivo PolyMag Reagent | Specifications

Package content	<p>IV-PN30500: 500µL of <i>In vivo</i> PolyMag</p> <p>IV-PN31000: 1mL of <i>In vivo</i> PolyMag</p> <p>IV-KC30210: 500µL of <i>In vivo</i> PolyMag + a Magnets set (IV-MAG1)</p> <p>IV-TK30210: 100µL of <i>In vivo</i> PolyMag + 1 cylinder magnet (ø 10mm)</p> <p>IV-MAG1, Magnet Set: 1 extra small cylinder (ø 2mm), 1 small cylinder (ø 5mm), 1 cylinder (ø 10mm), 1 square (18x18 mm) magnets</p> <p>IV-MAG2, Square Magnets set: 4 square magnets (18x18 mm)</p> <p>IV-MAG3, Cylinder Magnet set: 4 extra small cylinder (ø 2 mm), 4 small cylinder (ø 5 mm), 4 cylinder (ø 10 mm) magnets</p>
Shipping condition	Room Temperature
Storage condition	Store the <i>In vivo</i> PolyMag transfection reagent at +4°C upon reception
Shelf life	1 year from the date of purchase when properly stored and handled
Product Description	<i>In vivo</i> PolyMag is a cationic polymer-based magnetic nanoparticles formulation, designed for <i>in vivo</i> targeted transfection of nucleic acids.
Important notice	For research use only. Not for use in diagnostic procedures

Applications |

1. Nucleic acids

In vivo PolyMag has been developed for *in vivo* targeted transfection of various types of nucleic acids such as DNA, RNA. DNA/nanoparticles can be easily administrated through various injection routes such as systemic administration (intravenous, intra-artery) or local administration (intraperitoneal, intratumoral, intracerebroventricular, intramuscular).

Optimal conditions may vary depending on the nucleic acid, animal model, administration route and the target organ. Therefore, use the Table 1 as a starting point for DNA amount and volume of injection in mouse and rats.

Mouse			
Route of injection	Amount of nucleic acid	Total volume of injection according to animal weight	Site of injection
<i>Intravenous</i>	40 µg	200 µL (10-25 µL/g)	Tail vein
<i>Intramuscular</i>	10 to 100 µg	100 µL (50 µL x 2 sites of injection)	Caudal thigh
<i>Subcutaneous</i>	10 µg	200 µL (10-40 µL/g)	Scruff
<i>Intraperitoneal</i>	100 µg	400 µL (20 µL/g)	Lower Ventral Quadrants
<i>Intratumoral</i>	10 to 50 µg	100 µL (0.5 µL/mm ³)	Tumor
<i>Intracerebroventricular</i>	0.5 µg	2 µL	Brain ventricle
Rat			
Route of injection	Amount of nucleic acid	Total volume of injection	Site of injection
<i>Intravenous</i>	150 µg	2.5 mL (10-20 µL/g)	Tail vein, saphenous vein
<i>Intramuscular</i>	50 to 300 µg	300 µL (100 µL x 3 sites of injection)	Triceps, Quadriceps, Gluteals
<i>Subcutaneous</i>	5 to 10 µg	1.25 mL (5-10 µL/g)	Scruff, Back, Abdomen
<i>Intraperitoneal</i>	200 µg	2.5 mL (10-20 µL/g)	Lower Ventral Quadrants
<i>Intratumoral</i>	10 to 50 µg	100 µL (0.5 µL/mm ³)	Tumor
<i>Intracerebroventricular</i>	1 µg	10 µL	Brain ventricle

Table 1: Suggested amount of nucleic acid and volume of injection in mouse (20 g) and rat (250 g)

For more detailed protocols, see our Applications Notes on our website www.ozbiosciences.com or contact us at tech@ozbiosciences.com.

2. Magnets

Several kinds of magnets are provided with the *In vivo* PolyMag kit; use Table 2 to choose the best one adapted to your application.

Kind of magnet	Tissue
<u>Extra Small Cylinder</u> 2 mm (diameter) x 5 mm (height)	<ul style="list-style-type: none">• Brain area• Endothelial cells• Small tumors• Lymph node• Ovary• Adrenal gland
<u>Small Cylinder</u> 5 mm (diameter) x 5 mm (height)	<ul style="list-style-type: none">• Subcutaneous tumors• Salivary gland• Brain
<u>Cylinder</u> 10 mm x 5 mm (height)	<ul style="list-style-type: none">• Subcutaneous tumors• Pancreas• Spleen
<u>Square</u> 17 mm (length) x 17 mm (length) x 5mm (height)	<ul style="list-style-type: none">• Large organs• Large tumor• Muscle• Lung• Skin flap

Table 2: Examples of use of magnets

OZ Biosciences is currently proposing only those magnets. If you need specific magnet in terms of shape and size, please contact our technical service for obtaining fundamental properties of the magnet to purchase.

Protocol

Please refer to Table 1 to determine the required amount of DNA as well as volume injection. The DNA, *In vivo* PolyMag and injection solution should be at room temperature. We recommend using 1 μ L of *In vivo* PolyMag per μ g of DNA.

1. Reagent Preparation

- a. *In vivo* PolyMag. Before each use, vortex *In vivo* PolyMag vial. Add the required volume of *In vivo* PolyMag (according to DNA amount needed) to a sterile microtube.
- b. *DNA solution*. Dilute DNA in the final injection volume in a sterile vial (subtract the *In vivo* PolyMag volume).

2. Complexes formation.

- a. Add the DNA solution to the *In vivo* PolyMag and mix immediately by pipetting up and down.
- b. Incubate the complexes for 20 min at room temperature.

3. Injection.

- a. Place the magnet on your targeted tissue
- b. Slowly inject the complexes
- c. Let the magnet stand from 5 min to 1 h (Table 3 and next section).
Notes for intracerebroventricular or intra tumoral injections: Place the magnet few seconds after the complexes injection. Dye e.g. Fast Green FCF can be added to the complexes solution for a better monitoring of the injection.
- d. Monitor gene expression at the appropriate time point.

4. Magnetic incubation

The magnetic incubation time depends on the animal and the targeted tissue:

- for tumor, from 20 min (mouse, rat) to 1 hour (hamster, cat)
- for endothelial cells, from 5 to 20 min for mouse and rat, from 20 min to 1 h for rabbit or pig
- for peripheral tissue (e.g. stomach, gut, heart), 20 min
- for intracerebroventricular injection, 5 min

Refer to table 3, for other magnetic incubation times depending on target tissue, route of injection and magnet type.

Target tissue	Route of injection	Kind of magnet	Magnetic incubation
Tumor	Intravenous, Intra-arterial, Intratumoral	All kind	20 min to 1 h
Endothelial cells	Intravenous, Intra-arterial	Extra small Cylinder	5 min to 1 h
Heart	Intra-arterial	Cylinder	20 min
Liver	Intravenous	Cylinder, Square	10 min
Lung	Intravenous	Square	10 min
Pancreas	Intrapancreatic	Cylinder	20 min
Kidney	Intraperitoneal	Cylinder, Square	20 min
Gut	Ilea lumen	All kind	20 min
Stomach	Stomach lumen	Cylinder, Square	20 min
Brain	Intraventricular	Small Cylinder	5 min

Table 3: Suggested magnetic incubation time for various tissue

IMPORTANT NOTES:

- For long incubation time, (e.g. intratumoral injection), the magnet could be attached to the animal using adhesive tape in order to create a strong magnetic field in the area of the injection.
- Magnets can be easily handled with any magnetic surgical instruments (forceps, clamps, needle holders).
- Magnets can be sterilized by heat (steam sterilization or dry heat sterilization) or chemical means (ethanol 70%).

5. Bibliographic references

Please refer to the results sheet and to our website for a more exhaustive list of bibliographic references.

- Gupta A.K and Gupta M 2005 Synthesis and surface engineering of iron oxide nanoparticles for biomedical applications. *Biomaterials*. 26:3995-4021.
- Laurent N, Sapet C, Le Gourrierc L, Bertosio E and Zelphati O 2011 Nucleic acid delivery nanoparticles: the Magnetofection™ technology. *Therapeutic Delivery*. 2:471:482.
- Plank C, Zelphati O and Mykhaylyk O. 2011 Magnetically enhanced nucleic acid delivery. Ten years of magnetofection-progress and prospects. *Adv Drug Deliv Rev*. 63:1300-1331
- Ohashi K, Enomoto T, Joki Y, Shibata R, Ogura Y, Kataoka Y, Shimizu Y, Kambara T, Uemura Y, Yuasa D, Matsuo K, Hayakawa S, Hiramatsu-Ito M, Murohara T, Ouchi N. 2014 Neuron-derived neurotrophic factor functions as a novel modulator that enhances endothelial cell function and revascularization processes. *J Biol Chem*. 289:14132-44.

Related products for *in vivo* applications:

- **BrainFectIN** enables nucleic acids delivery into central nervous system of small animals.
- **In vivo DogtorMag** a cationic polymer-based magnetic nanoparticles formulation, designed for *in vivo* targeted transfection of nucleic acids.
- **In vivo SilenceMag** a cationic lipid-based magnetic nanoparticles formulation, designed to transfect small RNA, into target cell/ tissue *in vivo*.
- **In vivo ViroMag** an optimized nanoparticles formulation dedicated for *in vivo* transduction.

Purchaser Notification

Limited License

The purchase of the In vivo PolyMag grants the purchaser a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed this protocol). This reagent is intended for in-house research only by the buyer. Such use is limited to the transfection of nucleic acids as described in the product manual. In addition, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences. Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of the In vivo PolyMag. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, contact us at OZ Biosciences. Buyers may end this License at any time by returning all In vivo PolyMag reagents and documentation to OZ Biosciences, or by destroying all in vivo PolyMag components. Purchasers are advised to contact OZ Biosciences with the notification that a In vivo PolyMag is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s). This document covers entirely the terms of the In vivo PolyMag research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

In vivo PolyMag and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the kit components by following proper research laboratory practices.

EUROPE & ASIA OZ Biosciences SAS

163 avenue de Luminy
Case 922, zone entreprise
13288 Marseille cedex 09
France

Ph: +33 (0) 486 948 516
Fax: +33 (0) 463 740 015

contact@ozbiosciences.com
order@ozbiosciences.com
tech@ozbiosciences.com



USA & CANADA OZ Biosciences INC

7975 Dunbrook Road
Suite B
San Diego CA 92126
USA

Ph: + 1-858-246-7840
Fax: + 1-855-631-0626

contactUSA@ozbiosciences.com
orderUSA@ozbiosciences.com
techUSA@ozbiosciences.com

REV 06/2024