




Transfection reagent

Lullaby®

Lullaby siRNA Transfection Reagent


Protocol

Lullaby® Quick Protocol

To find your ideal silencing conditions with Lullaby®, we suggest to test increasing doses of siRNA (or miRNA): from 10 to 50nM per well.

Seed cells to be at 50-70% confluent the day of transfection

1



Prepare 3 tubes of siRNA (with different amounts of nucleic acids)*

2



96 well plate

24 well plate

6 well plate

10 nM/25nM/50nM in 50µL
serum-free medium or buffer*

10 nM/25nM/50nM in 50µL
serum-free medium or buffer*

10 nM/25nM/50nM in 100µL
serum-free medium or buffer*

Prepare 3 tubes of Lullaby® (with different amounts of reagent)*



3

96 well plate

24 well plate

6 well plate

0.5µL/1µL/1µL in 50µL of
serum-free medium or
buffer*

2µL/3µL/4µL in 50µL of
serum-free medium or
buffer*

8µL/10µL/14µL in 100µL of
serum-free medium or buffer*

Mix each tube of siRNA (step 2) to each tube of Lullaby® (step 3)*



4

96 well plate

24 well plate

6 well plate

siRNA *Lullaby®*

siRNA *Lullaby®*

siRNA *Lullaby®*

10nM + 0.5µL

10nM + 2µL

10nM + 8µL

25nM + 1µL

25nM + 3µL

25nM + 10µL

50nM + 1µL

50nM + 4µL

50nM + 14µL


Incubate 20 min at room temperature

5




Distribute each mix dropwise onto the cells

6




Incubate cells for 24 to 72h at 37°C until evaluation of transgene silencing

7



Choose the best ratio siRNA:Lullaby®

8



* Please refer to the following section "Important Notes"

IMPORTANT NOTES – Before you begin

- ✓ The siRNA optimal concentration required to achieve the best gene silencing effect depends highly on the cells, target and siRNA sequence; consequently, we suggest to first test a range of siRNA concentration from 10 to 50nM.
- ✓ For cell lines, seed the cells 24h before transfection in a 96-well plate, 24-well plate or 6-well plate in respectively 150 µL, 400 µL and 2 mL of complete culture medium.
- ✓ Allow reagents to reach RT and gently vortex them before forming complexes.
- ✓ For the preparation of Lullaby: it is important to add first the serum free medium to the tube and then add carefully the Lullaby reagent directly into the serum free medium without touching any plastic surface.
- ✓ Medium or buffer without serum & supplement must be used for the siRNA/Lullaby complexes preparation. Culture medium such as DMEM or OptiMEM or buffers such as HBS or PBS are recommended. In contrast, we do not recommend RPMI for preparing the complexes.
- ✓ We recommend respecting the order of addition of reagents; add the siRNA solution into the Lullaby solution.
- ✓ Dilute the reagent with deionized water for doses less than 1µL,

For additional information and protocols
(optimization, scaling, co-transfection...)
tips, troubleshooting or other applications



www.ozbiosciences.com

Any questions?



tech@ozbiosciences.com

Lullaby Reagent | Specifications

Package content	LL70500: 500µL of Lullaby LL71000: 1mL of Lullaby LL73000: 3 x 1mL of Lullaby
Shipping conditions	Room Temperature
Storage conditions	Store the Lullaby transfection reagent at +4°C upon reception
Shelf life	1 year from the date of purchase when properly stored and handled
Product description	Lullaby is the ideal siRNA transfection reagent for gene silencing
Important notice	For research use only. Not for use in diagnostic procedures

1. Cells Preparation

It is recommended to seed or plate the cells the day prior transfection, however cells can also be prepared few hours before the transfection or the reverse transfection. The suitable cell density will depend on the growth rate and the conditions of the cells. Best results are achieved if cells are at least 50-70 % confluent at the time of transfection (refer to Table 1).

Culture vessel	Number of adherent cells	Number of suspension cells	Cell overlay volume
96-well	$6 - 12 \times 10^3$	$4 - 8 \times 10^4$	100 μ L
24-well	$4 - 8 \times 10^4$	$25 - 50 \times 10^4$	400 μ L
6-well	$2 - 4 \times 10^5$	$1 - 2 \times 10^6$	1800 μ L

Table 1: Recommended number of cells to seed

2. siRNA/Lullaby complexes preparation

The siRNA and Lullaby solutions should have an ambient temperature and be gently vortexed prior to use.

- a. *siRNA solution*. Dilute the siRNA stock solution (for instance 1 μ M stock solution) in 50 or 100 μ L (refer to Table 2) of culture medium without serum and antibiotics.

Culture vessel	96-well		24-well		6-well	
Dilution serum-free medium	50μL		50 μL		100 μL	
Amount of siRNA (1 μM stock)*						
Final siRNA concentration	(μL)	(ng)	(μL)	(ng)	(μL)	(ng)
10 nM	2	27	5	67.5	20	270
20nM	4	54	10	135	40	540
50 nM	10	135	25	337.5	100	1350

* ng of siRNA was calculated on the basis of a MW = 13 500

Table 2: Suggested dilution procedure and amount of siRNA to test

- b. *Lullaby preparation*. Dilute 0.5 to 14 μ L of Lullaby in 50 or 100 μ L (refer to Table 3) of culture medium without serum and antibiotics.

Culture vessel	96-well	24-well	6-well
Dilution serum-free medium	50 μ L	50 μ L	100 μ L
Final transfection Volume	200 μ L	500 μ L	2 mL
Final siRNA concentration	Amount of Lullaby (μ L)		
10 nM	0.5	2	8
20nM	1	3	10
≥ 50 nM	1	4	14

Table 3: Recommended amount of Lullaby per nM of siRNA used

- c. Add the siRNA solution onto the Lullaby reagent and mix gently by carefully pipetting up and down.

- d. Incubate the mixture for 15-20 minutes at room temperature. Do not vortex or centrifuge!

3. Transfection

- a. Add the complexes drop by drop onto the cells and homogenize by gently rocking the plate side to side to ensure a uniform distribution of the mixture.
- b. Cultivate the cells at 37°C in a CO₂ incubator under standard conditions until evaluation of gene knockdown analysis.

NOTE: Depending on the siRNA amount, the gene targeted and the cell type, assays can be monitored 24 to 96h post-transfection.

Protocol | siRNA in suspension cells

1. Cell Preparation

The day before transfection split the cells at a density of 2 to 5 x 10⁵ cells / mL, so they are in excellent condition on the day of transfection. Incubate overnight in complete culture medium.

2. siRNA/Lullaby complexes preparation

The siRNA and Lullaby solutions should have an ambient temperature and be gently vortexed prior to use.

- a. *siRNA solution:* Dilute the siRNA stock solution (for instance 1 μM stock solution) in 50 or 100 μL (refer to Table 2) of culture medium without serum and antibiotics.
- b. *Lullaby preparation:* Dilute 0.5 to 14 μL of Lullaby in 50 or 100 μL (see Table 3) of culture medium without serum and antibiotics.
- c. Add the siRNA solution onto the Lullaby reagent and mix gently by carefully pipetting up and down.
- d. Incubate the mixture for 15-20 minutes at room temperature. Do not vortex or centrifuge!

3. Transfection

- a. While the complexes are incubating, prepare your cells in serum-free medium (or serum-containing medium) and transfer the appropriate volume to the culture dish according to Table 1. In 24-well plates for instance plate 2x10⁵ suspension cells just before transfection in 250 μL of serum free medium. Generally, serum-free condition leads to higher transfection efficiency.
- b. Next, add the complexes directly onto the cells dropwise and all over the well. Important: gently mix complexes with the cells by pipetting the culture medium up and down (3-4 times to disrupt potential cell clumps and to ensure contact of the complexes with cells
- c. Incubate 3 to 6 h (4h is commonly used) in serum-free medium at 37°C under 5% CO₂.
- d. If transfections are performed in serum free medium, add serum to adjust its concentration.
- e. Cultivate the cells at 37°C in a CO₂ incubator under standard conditions until evaluation of gene knockdown analysis.

NOTE: Depending on the siRNA amount, the gene targeted and the cell type, assays can be monitored 24 to 96h post-transfection

In order to get the best out of **Lullaby® siRNA transfection reagent**, several parameters can be optimized:

- Ratio of Lullaby reagent to siRNA
 - siRNA dose used, which strongly depends on the efficiency and specificity of your siRNA
 - Cell type, cell density and incubation time
1. Start by optimizing the ratio Lullaby / siRNA. To this end, use a fixed amount of siRNA and vary the amount of **Lullaby®** as detailed in the Table 4. The reagents can be pre-diluted in deionized water and aliquots of the resulting dilutions are incubated with siRNA. Diluted Lullaby solution has to be freshly prepared.

Culture vessel	96-well	24-well	6-well
Dilution serum-free medium	50 µL	50 µL	100 µL
Final transfection Volume	200 µL	500 µL	2 mL
Final siRNA concentration	Amount of Lullaby (µl)		
5 nM	0.25 – 0.5 – 1 – 1.5	0.5 – 1 – 2 - 3	2 – 4 – 6 - 8
10 nM	0.25 – 0.5 – 1 – 1.5	1 – 2 – 3 - 4	4 – 8 – 12 - 16
20nM	0.5 – 1 – 2 - 3	1.5 – 3 – 4 - 6	7 – 10 – 15 - 20
≥ 50 nM	0.5 – 1 – 2 - 3	2 – 4 – 6 - 8	10 – 14 – 18 - 22

Table 4: Recommended amount of Lullaby per nM of siRNA used

2. Thereafter, optimize the siRNA dose with the fixed ratio of Lullaby / siRNA that has been previously optimized (refer to Table 4).
3. After having identified the optimal quantity of Lullaby reagent and siRNA, you could pursue the process by optimizing the cell number (density) and time course of your experiment.

Additional products for your silencing experiments

- **Lullaby Stem** for siRNA transfection into stem cells
- **SilenceMag** Magnetofection reagent dedicated to siRNA transfection into hard-to-transfect cells

Purchaser Notification

Limited License

The purchase of the Lullaby kit grants the purchaser a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in this protocol). This reagent is intended for in-house research only by the buyer. Such use is limited to the transfection of nucleic acids as described in the product manual. In addition, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences. Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of the Lullaby kit. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, contact us at OZ Biosciences. Buyers may end this License at any time by returning all Lullaby kit reagents and documentation to OZ Biosciences, or by destroying all Lullaby components. Purchasers are advised to contact OZ Biosciences with the notification that a Lullaby kit is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s). This document covers entirely the terms of the Lullaby kit research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

Lullaby kit and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the kit components by following proper research laboratory practices.

EUROPE & ASIA OZ Biosciences SAS

163 avenue de Luminy
Case 922, zone entreprise
13288 Marseille cedex 09
France

Ph: +33 (0) 486 948 516
Fax: +33 (0) 463 740 015

contact@ozbiosciences.com
order@ozbiosciences.com
tech@ozbiosciences.com



USA & CANADA OZ Biosciences INC

7975 Dunbrook Road
Suite B
San Diego CA 92126
USA

Ph: + 1-858-246-7840
Fax: + 1-855-631-0626

contactUSA@ozbiosciences.com
orderUSA@ozbiosciences.com
techUSA@ozbiosciences.com