

## Description

Ready-to-use stabilized **DiL-labeled LNPs** for *in vitro* cellular and subcellular distribution and *in vivo* imaging.

**Concentration:** LNPs containing 1 mole % of DiL

**Buffer:** PBS, 10 % sucrose

LNPs represent one of the most effective and safe delivery systems for the translational success of nucleic acid drugs. LNPs are lipidic spherical vesicles formed by a combination of four main components: an ionizable cationic lipid, a helper phospholipid, cholesterol & a pegylated lipid, each having distinct functions. LNPs not only protect RNA from degradation, but also facilitate intracellular uptake and thus potentiate its efficacy.

**NanOZ Fluo LNP-DiL (Fig.1)** is designed for the effective monitoring and imaging of intracellular trafficking and biodistribution of LNPs *in vitro* and *in vivo* respectively.

**NanOZ Fluo LNP-DiL** contains a lipophilic carbocyanine dye DiL having an excitation max. at 549 nm, and emission max. at 565 nm.




Fig.1. Schematic representation of LNPs-DiL

## Applications

### DiL detection

The DiL has an excitation peak at 549 nm and emission peak at 565 nm and can be monitored by flow cytometer, plate readers and confocal microscopy.

## Quality Controls

Items	Specification	Standard QC	Superior Grade QC*
Identity	Size	✓	✓
	Charge	✓	✓
Safety	Sterility	✓	✓
	Endotoxin		✓
	Mycoplasma detection		✓
Characterization	Lipid content		✓

\* Contact us to get a quote.

Certificate of analysis on demand.

## Use, handling and storage

*For Research Use Only. Not for use in humans. Not for use in diagnostic or therapeutic purposes.*

**Long term storage (6 months):** -80°C.

**Short term storage (2 months):** +4°

We recommend minimizing freeze-thaw cycles to preserve LNPs integrity.

## Kit contents

**LFL11000:** 1mL (4\*250µL) of LNP-DiL, (4.5 mM) having composition as described in table below:

Lipid mix components	Molecular weight	Molar ratio
FP105	731.2	50.0
DSPC	790.2	9.0
Cholesterol	386.7	38.5
DMG-PEG 2000	2509.2	1.5
DiL	933.9	1.0
Total		100.0

## Related Products

Ref	Description
#LNP10500	Empty LNP
#LNP10500mRNA11	NanOZ-LNP/mRNA (GFP).
#LNP10500mRNA12	NanOZ-LNP/mRNA (LUC).
#LNP10500mRNA41	NanOZ-LNP/mRNA (OVA).

Custom LNPs & mRNAs are also available now!

### Purchaser Notification | Conditions of Sale

This product is sold in accordance with our general conditions of sale that you can find on our website: <https://ozbiosciences.com/content/3-terms-and-conditions>.