

To find the ideal conditions, HYPE-293™ SG reagent can be tested at ratios **1:1**, **2:1** and **3:1** (respectively 1 µL, 2 µL and 3 µL of HYPE-293 SG per µg of DNA). DNA quantities vary from **0.75 µg** to **1.5 µg** of DNA per mL of culture medium. B293 SG reagent can be used at 1/100 or 1/200 of the total volume.


1. Prepare the cells



2. Prepare 3 identical tubes of DNA



3. Prepare 3 tubes of HYPE-293 SG corresponding to 3 different ratios



4. Mix each DNA suspension to each tube of HYPE-293 SG dilution



5. Transfect cells and add B293 SG reagent



6. Incubate cells under orbital shaking at 37°C until evaluation of protein production

Chose the optimal conditions (DNA quantity, HYPE-293 SG ratio, B293 SG volume).

Before you begin

HYPE-293 Superior Grade (SG) is compatible with most culture media for protein production except for CD293 from Life Technologies. Do not use culture medium containing high antibiotic level (up to 0.5 X penicillin/streptomycin final concentration) or high Pluronic® surfactant concentration (up to 0.01% w/v final concentration) to avoid dramatic impact on protein production level.

The use of B293 SG reagent is highly recommended yet optional.

HYPE-293 Superior Grade is identical in synthesis and formulation to HYPE-293 Research Use Only Grade (RUO) but provides more quality controls, batch validation and documentation as required for preclinical trials.

Protocol | DNA Transfection

1. Cells preparation:

Cell culture maintenance: sub-culture the cells at a density of 0.5-2x10⁶ cells/mL for each passage (48-72 h). Avoid high cell density and keep cell growth conditions consistent for reproducibility.

18-24 h before transfection, dilute the cells to 0.6-0.8x10⁶ cells/mL and incubate on orbital shaker (~125 rpm) at 37 °C, 8% CO₂.

The day of transfection, dilute the cells to 1x10⁶ cells/mL (cell density should be about 1.2-1.5x10⁶ cells/mL). Transfer the volume of cells needed as described in Table 1.

Cell Culture			DNA		HYPE-293 SG		B293 SG
10 ⁶ cells per mL			1.5 µg/mL		2 µl per µg DNA		1X final
Vol.	Flask.	Cell nb.*	µg	Vol.	µL	Vol.	µL
1 mL	NA	1x10 ⁶	1.5 µg	50 µL	3 µL	50 µL	11 µL
30 mL	125 mL	30x10 ⁶	45 µg	0.5 mL	90 µL	0.5 mL	310 µL
250 mL	1 L	250x10 ⁶	375 µg	5 mL	750 µL	5 mL	2.6 mL
1 L	3 L	1x10 ⁷	1.5 mg	20 mL	3 mL	20 mL	10.4 mL

* The day of transfection cell density should be at 1 x 10⁶ cells/mL.

Table 1: Suggested volumes of HYPE-293 SG, B293 SG and DNA quantity

2. HYPE-293 SG/DNA complexes preparation:

a. **HYPE-293 SG:** vortex the reagent and dilute the indicated quantity of HYPE-293 SG in 50 µL to 20 mL of culture medium without serum and supplement.

b. **DNA:** dilute the indicated quantity of DNA in 50 µL to 20 mL of culture medium w/o serum and supplement.

c. **Complexes:** add the DNA solution to the HYPE-293 SG solution and mix gently by carefully pipetting up and down.

Incubate the mix at room temperature for 20 minutes.
Do not vortex or centrifuge.

Refer to table 1 for volumes depending on size of cell culture dishes.

3. Transfection

a. Add the HYPE-293 SG/DNA complexes dropwise into cell culture bottle while gently swirling the flask to ensure a uniform distribution.

b. Add the B293 SG reagent – 1X final directly into the vessel containing cells.

c. Cultivate the cells under standard conditions (~125 rpm, 37°C, 8% CO₂) for 1 to 7 days depending on the type of protein expression. No medium change is required during the incubation period.

Our technical team is at your disposal for any questions or optimization procedures:

✉ - tech@ozbiosciences.com

Use, handling and storage

For Research Use Only. Not for use in humans. Not for use in diagnostic or therapeutic purposes.

Shipping conditions: Room Temperature

Storage conditions HYPE-293 SG & B293 SG: -20°C

Shelf life: 1 year from the date of purchase

We recommend minimizing freeze-thaw cycles to preserve HYPE-293 SG integrity.

Quality control:

	Hype 293	Hype 293 SG
Sterility	✓	✓
Biological activity	✓	✓
Chemical characterization		✓
Purity		✓
Absence of Mycoplasma*		✓
Endotoxin ≤ 1 EU/mL**		✓

*USP <63> Mycoplasma – **USP <85> Bacterial Endotoxin test / Certificate of analysis on demand.

Kit contents

HY293SG15: 1.5 mL of HYPE-293 + 5 mL of B293

HY293SG30: 2 x 1.5 mL of HYPE-293 + 2 x 5 mL of B293

HY293SG150: 15 mL of HYPE-293 + 50 mL of B293

HY293SG300: 2 x 15 mL of HYPE-293 + 2 x 50 mL of B293

Related Products

Ref	Description
#HY29315	HYPE-293, bioproduction in HEK-293-S

Purchaser Notification | Conditions of Sale

This product is sold in accordance with our general conditions of sale that you can find on our website: <https://ozbiosciences.com/content/3-terms-and-conditions>.