

Content and product information

Formulation: AlPO₄, Aluminum phosphate gel

Appearance: White gelatinous precipitate

AlumVax Phosphate is wet gel (colloidal) of aluminum phosphate 2%, provided as a ready-to-use suspension. It is sterilized and aseptically filled. AlumVax phosphate is an amorphous aluminum hydroxyphosphate which is negatively charged at physiological pH (pI=5–7), suitable for adsorption of positively charged or neutral, alkaline proteins.

Description

Aluminum salts are the most common adjuvant used in approved prophylactic vaccines because of their excellent safety profile and ability to enhance protective humoral immune response. Since more than 80 years, it has been observed that aluminum compounds act by a depot effect and also by direct activation of the immune cells.

Adsorption or entrapment of antigens in aggregates through hydrophobic and electrostatic interactions favors a high local antigen concentration and improved uptake by antigen presenting cells (APC).

Alum Phosphate stimulates Th2 response through the release of Th2-associated cytokines (IL4, IL-5, IL-13...) and Th2-associated antibodies (IgG1 & IgE). It increases Ag-specific CD4+ T Cell proliferation and promotes NALP3 inflammasome activation and caspase 1-mediated release of IL-1 and IL-18.

NOTE: Alum is frequently used as an alternative to Freund's adjuvants, as it is less hazardous and less likely to cause tissue necrosis at the injection site.

Use, handling and storage

For Research Use Only. Not for use in humans. Not for use in diagnostic or therapeutic purposes.

Shipping conditions: Room Temperature

Storage conditions: Room Temperature

Shelf life: 1 year from the date of purchase

⚠ Do not freeze.

Kit contents

AP0050: 50 mL of AlumVax Phosphate 2%

AP0250: 250 mL of AlumVax Phosphate 2%

Certificate of analysis on demand.

Method | Protocol

Recommendations before starting:

The inoculum should be free of extraneous microbial contamination; filtration of the antigen before mixing with the adjuvant is recommended.

Preparation of AlumVax Phosphate-immunogens mix is easier than other adjuvants as it does not require laborious emulsification.

1. Ensure complete re-suspension of AlumVax phosphate adjuvant by vigorous shaking of the vial before use.

2. Dilute antigen mixture in saline buffer or phosphate buffer for a final immunogen concentration of **10-100 µg/100 µL**.

3. Mix AlumVax phosphate adjuvant with an equal volume of antigen solution for a **1:1** ratio:

a. Add AlumVax phosphate dropwise with constant mixing to the immunogen solution.

b. Pipet up and down several times to ensure correct absorption of antigen by alum adjuvant and incubate **5-10 min**.

NOTE: Ratio can be optimized from **1:1** (100 µL adjuvant per 100 µL antigen) to **1:9** (100 µL adjuvant per 900 µL antigen).

4. Inject into the animal according to table 1 below; the volume depends on the site of injection.

Typical routes of administration include subcutaneous (SC), intramuscular (IM), intradermal (ID) or intraperitoneal (IP).

Species	Max vol. / Site	Primary Injection	Subsequent Injection(s)
Mice, hamsters	0.1 mL	SC	SC
Mice, hamsters	0.05 mL	IM ^Δ	IM ^Δ
Mice	0.5 mL	IP [×]	SC, IM ^Δ
Guinea pigs, rats	0.2 mL	SC, IM ^Δ	SC, IM ^Δ
Rabbits	0.25 mL	SC, IM	SC, IM
Rabbits	0.025 mL	ID	SC, IM
Sheep, goats, donkeys, pigs, chickens	0.5 mL	SC, IM	

Δ Not recommended in general, particularly for viscous adjuvants.

× Not recommended for pAb production.

Table 1: Maximum volumes for injection of immunogen/adjuvant mixtures per site of injection for different animal species (Adapted from Leenars MPPA, Hendriksen CFM et al., 19).

Purchaser Notification | Conditions of Sale

This product is sold in accordance with our general conditions of sale that you can find on our website: <https://ozbiosciences.com/content/3-terms-and-conditions>.

Results

Results presented below demonstrate the effect of Aluminum Phosphate adjuvant on immune system response:




Figure 1. Adjuvant effect of AIPO4 for HBV DNA vaccines. Aluminum phosphate had a powerful adjuvant effect for Hepatitis B (HBV) DNA vaccines in mice (adapted from Wang S. *et al.*, *Vaccines*. 2000; 18:1227-35).




Figure 2. Efficient delivery of protein antigen. BALB/C mice were vaccinated i.m. with recombinant Hepatitis B surface antigen with or without AIPO4. 4 weeks sera were tested for anti-HBsAg antibodies by ELISA (adapted from Kwissa M. *et al.*, *J Mol Med*. 2003; 81:502-510).

Related Products

Ref	Description
#AH0250	AlumVax Hydroxide 2%
#SQ0100	SqualVax, squalene oil-in-water emulsion
#IFA0100	Incomplete Freund's Adjuvant (IFA)
#CFA0100	Complete Freund's Adjuvant (CFA)
#CV02000	CaLiVax-DOTAP Adjuvant
#LV02000	LipoVax NTA(Ni)