



Viral Applications

ViroMICST Stem

Transduction enhancer

Stem cells infection during magnetic cell purification


Protocol

Quick proto schema for infection/transduction on separation column

The following protocol is given for MS column with a void volume of 60 μL ; adjust or adapt volumes depending on the column format. If the lentiviral transduction/infection conditions are unknown, we recommend starting with a MOI of 2 using a lentiviral vector encoding for a fluorescent protein and 3 different volumes of ViroMICST Stem Reagent ST with 3 μL of Reagent VM.


ViroMICST STEM Working Solution

Before Beginning, prepare 3 ViroMICST Stem Working Solutions (WS) by mixing 2, 5 or 10 μL Reagent ST with 3 μL Reagent VM



	ViroMICST Stem working solution	
	Reagent ST	Reagent VM
Tube 1	-	-
Tube 2	2 μL	3 μL
Tube 3	5 μL	3 μL
Tube 4	10 μL	3 μL


Magnetically label cells, refer to instruction manual




1a MAGNETICALLY LABELLING OF TARGET CELL POPULATION

Incubation time: 25 min.


1b COMPLEXES FORMATION




Refer to instruction manual



2 LOAD CELLS ON COLUMN



3 LOAD COMPLEXES ON COLUMN
Incubation time: 30 min.



4 FLUSH THE CELLS


5 ASSAY (24/96H)

Quick proto schema for infection/transduction on column-free separation magnet

The following protocol is given for cell sorting using column-free immunomagnetic separation magnet in 2.5 - 5 mL volume; adjust or adapt volumes depending on cell sorting apparatus format or cell culture volume. If the lentiviral transduction/infection conditions are unknown, we recommend starting with a MOI of 2 using a lentiviral vector encoding for a fluorescent protein and 3 different volumes of ViroMICST Stem Reagent ST with 3 μ L of Reagent VM


ViroMICST STEM Working Solution

Before Beginning, prepare 3 ViroMICST Stem Working Solutions (WS) by mixing 2, 5 or 10 μ L Reagent ST with 3 μ L Reagent VM




	ViroMICST Stem working solution	
	Reagent ST	Reagent VM
Tube 1	-	-
Tube 2	2 μ L	3 μ L
Tube 3	5 μ L	3 μ L
Tube 4	10 μ L	3 μ L


Magnetically label cells, refer to instruction manual




1a MAGNETICALLY LABELLING OF TARGET CELL POPULATION




Incubation time: 25 min.




2 INCUBATE LABELED CELLS IN MAGNET



3 POUR OFF UNDESIRED CELLS AND RESUSPEND CELLS IN SAME INITIAL VOLUME



4 ADD MAGNETIC COMPLEXES Incubation time: 30 min.




5 REMOVE TUBE FROM MAGNET & INCUBATE CELLS FOR 24-96H

IMPORTANT NOTES – Before you begin


- ✓ This short protocol is suitable for transducing 0.1 to 1×10^7 labeled-cells per immunomagnetic cell separation (magnet or column)
- ✓ Viro-MICST Stem does not interfere with magnetic cell sorting devices and we recommend following rigorously the cell sorting protocol given by the manufacturer.
- ✓ Using column purification, high yield of purification is reached when using at least two columns; apply Viro-MICST always on the last one of the process.
- ✓ Allow reagents to reach room temperature before starting
- ✓ Do not change your protocol for cell sorting after immunomagnetic labelling, simply add complexes of viral particles/ViroMICST Stem to the cell sorting apparatus (column or magnet).
- ✓ Do not mix Reagent ST or Reagent VM with any other viral enhancer or chemicals
- ✓ If transduction conditions are unknown, we recommend beginning with a MOI of 2
- ✓ We suggest starting with 2, 5 and 10 μL of ViroMICST STem Reagent ST with 3 μL Reagent VM – on the basis of MS column size from Miltenyi or 2.5 to 5 mL volumes for cell sorting magnet
- ✓ Refer to optimization protocol for gaining more information on how to optimize conditions.

OZ Biosciences does not provide magnetic cell separation systems ; OZ Biosciences does not provide the immunomagnetic beads. We recommend following rigorously the cell sorting protocols for both items given by the manufacturers.

For additional information and protocols (optimization, scaling, co-transfection...) tips, troubleshooting or other applications

 www.ozbiosciences.com

Any questions?

 tech@ozbiosciences.com

ViroMICST Stem Transduction Enhancer | Specifications

Package content	VMXS100: ViroMICST Stem 100 transductions (500 µL of Reagent ST + 300 µL of Reagent VM) VMXS300: ViroMICST Stem 300 transductions (1500 µL of Reagent ST + 900 µL of Reagent VM)
Shipping condition	Room Temperature
Storage conditions	Upon reception and for long-term use, store Reagent ST at -20°C and Reagent VM at +4°C.
Shelf life	1 year from the date of purchase when properly stored and handled
Product description	Viro-MICST Stem Transduction Enhancer was specifically designed for stem cell lentiviral transduction in order to achieve high infection rate directly on magnetic cells sorting devices (column or magnet)
Important notice	For research use only. Not for use in diagnostic procedures

1. Cells Preparation

Refer to columns manufacturer for target cells pre-enrichment detailed protocol as well as for labelling cells with nano-sized or micro-sized immunomagnetic beads.

2. Optimization Protocol

As a starting point, we recommend trying 3 different volumes of Reagent ST with 3 μ L of Reagent VM ; however several parameters can be optimized such as : MOI, volumes of each tube composing ViroMICST Stem per column, cell number...

It is generally suggested to optimize one parameter at a time while keeping the other parameters constant.

- a. Optimizing ViroMICST Stem volume – example for 0.1 to 1×10^7 labeled-cells per immunomagnetic cell separation (magnet or column) : using a fixed MOI, vary the volume of each tubes composing ViroMICST Stem working solution.

NOTE : These conditions are given for MS column or using a separation magnet for up to 5 mL cell suspension. Adapt volumes according to other cell sorting columns format or separation volumes.

- b. Quantity of viral particles: once the ideal ratio of ViroMICST Stem is found, keep it fixed and vary MOI.
- c. Cell number: finally use the best ratio and MOI found in (a) and (b) and vary the cell number to be assayed.

NOTES

Additional products for Viral Transduction Enhancement

- **ViroMag** for enhancing viral transduction efficiency (suitable for all viruses)
- **ViroMag R/L** for enhancing Lentiviral and Retroviral transduction efficiency
- **AdenoMag** specific for Adenoviral and AAV transduction

Additional products for Virus Production

- **Helix-IN DNA Transfection Reagent**

Additional products for Virus Capture and Concentration

- **Mag4C-LV** for Lentiviruses
- **Mag4C-AD** for Adenoviruses

Purchaser Notification

Limited License

The purchase of the ViroMICST Stem reagent grants the purchaser a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in section 1, Kit Contents) for the sole purpose of in-house research only, provided that no license, right or permission is granted hereunder to a non-academic, for-profit or commercial Licensee to use the ViroMICST Stem reagent for ex vivo gene therapy for hemoglobinopathies. The license does not include the use for any commercial or development purpose, including but not limited to any use for a) manufacturing, production, quality control, b) providing services, information or data, c) therapeutic, diagnostic, vaccine or prophylactic purposes or d) any applications which require regulatory approval as well as e) any clinical activities *in vivo* or *ex vivo*. The licensed use is limited to transfection of nucleic acids as described in the product manual. In addition, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences.

To inquire about separate licenses for the express purpose of non-research use or applications of the ViroMICST Stem reagent, or to obtain authorization to transfer or use the enclosed material, please contact the Director of Business Development at OZ Biosciences. Buyers may end this License at any time by returning all ViroMICST Stem reagent material and documentation to OZ Biosciences, or by destroying all ViroMICST Stem reagent components. Purchasers are advised to contact OZ Biosciences with the notification that a ViroMICST Stem reagent kit is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s). This document covers entirely the terms of ViroMICST Stem reagent research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

The ViroMICST Stem reagent and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the kit components by following proper research laboratory practices.

EUROPE & ASIA OZ Biosciences SAS

163 avenue de Luminy
Case 922, zone entreprise
13288 Marseille cedex 09
France

Ph: +33 (0) 486 948 516
Fax: +33 (0) 463 740 015

contact@ozbiosciences.com
order@ozbiosciences.com
tech@ozbiosciences.com



USA & CANADA OZ Biosciences INC

7975 Dunbrook Road
Suite B
San Diego CA 92126
USA

Ph: + 1-858-246-7840
Fax: + 1-855-631-0626

contactUSA@ozbiosciences.com
orderUSA@ozbiosciences.com
techUSA@ozbiosciences.com

Rev. 03/2021 - CT