



# Custom mRNA Service

INTEGRATED MRNA SERVICE FROM GENE SYNTHESIS  
TO IVT MRNA PRODUCTION





# CUSTOM mRNA SERVICE

With over 15 years of nucleic acid delivery expertise, OZ Biosciences is providing you with the **best-in-class** service for the production of **high quality mRNAs**.



A very flexible platform with a lot of options you can choose from:



We ship your **fully custom-made mRNA**

Don't hesitate to ask for other options such as **concentration, buffer, desired UTR, codon optimization** etc...

# READY-TO-USE mRNA

Category	Product name
<b>Reporter gene mRNA</b>	β-Gal mRNA
	mCherry mRNA
	GFP mRNA
	Tomato mRNA
	F-Luc mRNA
	R-Luc mRNA
<b>Fluorescent mRNA</b>	Cy5 F-Luc mRNA
	Cy5 GFP mRNA
<b>Gene editing mRNA</b>	Cas9 mRNA
	CRE mRNA
	Cas13d mRNA
<b>Vaccine / antigen mRNA</b>	OVA mRNA
	Spike DELTA mRNA
	Spike Omicron mRNA
	N SARS Cov2 mRNA
	HA-H1N1 mRNA
	HA-H2N3 mRNA
	Spike SARS CoV-2 mRNA
<b>Gene replacement mRNA</b>	EPO mRNA
<b>Interleukin mRNA</b>	h IL-2 mRNA

**OZ Biosciences** offers mRNAs that mimic fully processed mature mRNAs.

These mRNAs are stabilized with **5' Cap 1** structure and **3' poly(A) tail** and are optimized to yield improved stability & performance.

They are modified with **5-methoxyuridine** (5moU replaces U) or **N1-methyl-pseudouridine** to reduce innate immune responses. We also offer **unmodified mRNAs**.



**OZ BIOSCIENCES**  
The art of delivery systems

## USA & CANADA

OZ Biosciences USA Inc.  
7975 Dunbrook Road, Suite B,  
San Diego, CA 92126  
USA  
Ph: +1 858 246 7840  
Fax: +1 855 631 0626  
orderusa@ozbiosciences.com

## FRANCE & EUROPE

OZ Biosciences SAS  
Parc Scientifique de Luminy  
zone entreprise, case 922  
13288 Marseille Cedex 09  
FRANCE  
Ph: +33 486 94 85 16  
Fax: +33 486 94 85 15  
order@ozbiosciences.com