



**MATERIAL SAFETY DATA SHEET**

**SECTION 1 - PRODUCT AND COMPANY INFORMATION**

• **Product name: NanOZ LNP-DIY-MULTI-ORGAN**

- **Catalog number: LDIY022.**
- **Chemical name or synonyms: None**
- **Contact:**

OZ Biosciences  
 Parc Scientifique de Luminy – Zone Luminy Entreprises  
 163 Avenue de Luminy – Case 922  
 13288 Marseille cedex 9, France  
 Phone: +33 (0)4.86.94.85.16 / Fax: +33 (0)4.86.94.85.15  
[www.ozbiosciences.com](http://www.ozbiosciences.com) / [contact@ozbiosciences.com](mailto:contact@ozbiosciences.com)

**SECTION 2 – HAZARDS IDENTIFICATION**

**GHS Classification**

Not classified as Hazardous

**GHS Label element**

Signal Word: None  
 Hazard Statements: None  
 Pictograms: None

This product contains hazardous ingredients (Yes; No; Not Applicable or no information): **No.**

This product contains hazardous carcinogens (Yes; No; N/A): **No**

This product contains hazardous ingredients to the environment. (Yes; No; N/A): **No**

**SECTION 3 – COMPOSITION, INFORMATION ON INGREDIENTS**

**Description:** Dried film made from a proprietary composition of lipid nanoparticles (LNPs).

Chemical	EC No	CAS No	Weight-%	Classification according to Regulation (EC) No. 1272/2008 [CLP]	REACH RegistrationNumber
OZB2 Ionizable lipid	Not Listed	Not Assigned	<2	Not Listed	
Cholesterol	200-353-2	57-88-5	<1	Not Listed	
1,2-Dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (DMG-PEG)	Not listed	160743-62-4	<1	Not Listed	
1,2-Distearoyl-sn-glycero-3-phosphocholine	212-440-2	816-94-4	<1	Not Listed	

Ingredients are not hazardous or their concentrations do not exceed the limit specified according to 29 CFR 1910.1200 Hazard Communication Standard (USA) and Directive 1999/45/EC-2001/59/EC (EU).

\*\*This is proprietary and confidential information.

**SECTION 4 – FIRST AID MEASURES**

**Effects and symptoms:**

- *Inhalation* May be harmful if inhaled. Material may be irritating to mucous membranes and upper respiratory tract.
- *Ingestion* May be harmful in case of ingestion.
- *Skin contact* May cause skin irritation, redness, itching.
- *Eye contact* Slightly hazardous in case of eye contact.
- *Aggravating conditions* Chronic exposure may develop allergic reactions in sensitive individuals.

**First-Aid measures:**

- *Inhalation* If inhaled, remove to fresh air. If breathing is difficult, give oxygen. If not breathing, give artificial respiration. Get medical attention immediately.
- *Ingestion* If swallowed, wash out mouth with water provided person is conscious. Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Get medical attention if symptoms appear.
- *Skin contact* In case of contact, immediately wash skin with soap and copious amounts of water.
- *Eye contact* In case of contact, immediately flush eyes with plenty of water for at least 15 minutes. Get medical attention immediately.
- *Notes to physician* Not available.
- *Protection of first-aiders* Not available.

**SECTION 5 – FIRE FIGHTING MEASURES**

- **Extinguishing media:** Water spray, carbon dioxide, dry chemical powder or appropriate foam.
- **Flammability of the product:** Not available.
- **Flash Point:** Not available.
- **Fire hazards in presence of various substances:** Not considered to be flammable.
- **Fire fighting media and instructions:** Use an extinguish agent suitable for surrounding fires.
- **Protective clothing (fire):** Be sure to use an approved/certified respirator or equivalent.
- **Hazardous thermal decomposition products:** Not available.

**SECTION 6 – ACCIDENTAL RELEASE MEASURES**

- **Personal precautions:** Safety glasses, lab coat and gloves.
- **Environmental precautions and clean-up methods:** Absorb with an inert material and put the spilled material in an appropriate waste disposal. Finish cleaning by spreading water on the contaminated surface and allow evacuating through the sanitary system.
- **Small spill and leak:** Dilute with water and mop up, or absorb with an inert dry material and place in an appropriate waste disposal container. Finish cleaning by spreading water on the contaminated surface and dispose of according to local and regional authority requirements.
- **Additional information:** No dangerous substances are released.

**SECTION 7 – HANDLING AND STORAGE**

- **Handling:** Avoid contact with skin, eyes and clothing. Avoid prolonged or repeated exposure.
- **Storage:** Keep container tightly closed. Keep container in a cool, well-ventilated area. Store at - 20°C for long-term storage.
- **Intended use:** Refer to the instruction booklet for proper and intended use. Otherwise, contact supplier for specific application.
- **Packaging materials suitable:** Use original container.

**SECTION 8 – EXPOSURE CONTROLS, PERSONAL PROTECTION**

- **Engineering Controls:** Use process enclosures, local exhaust ventilation, or other engineering controls to keep airborne levels below recommended exposure limits. Ensure that eyewash stations and safety showers are proximal to the work-station location.
- **Personal protection**
  - **Eyes:** Safety glasses.
  - **Body:** Lab coat and gloves.
  - **Respiratory:** Respirator is not needed under normal and intended conditions of use, if exposures are kept below established limits.

**SECTION 9 – PHYSICAL AND CHEMICAL PROPERTIES**

- Appearance
  - **Physical State:** Solid film
  - **Color:** Colorless
  - **Odor:** Not available.
  - **Taste:** Not available.
- **Molecular weight:** Not available.
- **pH:** Not Available
- **Boiling point:** Not Available
- **Decomposition temperature:** Not available
- **Melting point:** Not available.
- **Flash point:** Not available.
- **Dispersion properties:** Not applicable.
- **Solubility:** Readily dispersable in water.

## SECTION 10 – STABILITY AND REACTIVITY

- **Stability and reactivity:** The product is stable under recommended storage conditions.
- **Conditions to avoid:** Exposure to light, air or moisture over prolonged periods.
- **Materials to avoid:** Strong acids, strong bases, strong oxidizing agents.
- **Hazardous polymerization:** Will not occur.
- **Hazardous decomposition products:** Will not occur.

## SECTION 11 – TOXICOLOGICAL INFORMATION

## 11.1. Information on toxicological effects

- **General Information:** Toxicological properties have not been thoroughly investigated. The following information is available for the individual ingredients.
- **Known Clinical Effects:** Based on clinical trials in humans, possible adverse effects following intravenous exposure to this compound may include: muscle pain, abnormal redness of skin (erythema), fever, and sleep disturbances.
- **Acute Toxicity: (Species, Route, End Point, Dose)**  
N/A
- **Irritation / Sensitization: (Study Type, Species, Severity)**  
N/A
- **Reproduction & Development Toxicity: (Duration, Species, Route, Dose, End Point, Effect(s))**  
N/A
- **Genetic Toxicity: (Study Type, Cell Type/Organism, Result)**  
N/A
- **Carcinogenicity See below**

**Cholesterol**

IARC Group 3 (Not Classifiable)

## SECTION 12 – ECOLOGICAL INFORMATION

**Environmental Overview:** Environmental properties have not been investigated. Releases to the environment should be avoided.

## 12.1. Toxicity Aquatic Toxicity: (Species, Method, End Point, Duration, Result)

N/A

## 12.2. Persistence and degradability

Persistence and degradability No information available.

## 12.3. Bioaccumulative potential

Bioaccumulation No information available.

## 12.4. Mobility in soil

Mobility in soil No information available.

## 12.5. Results of PBT and vPvB assessment

PBT and vPvB assessment No information available.

## 12.6. Other adverse effects

Other adverse effects No information available.

## SECTION 13 – DISPOSAL CONSIDERATIONS

- **Waste stream:** Not available.
- **Waste disposal method:** Contact a licensed professional waste disposal service to dispose of this material. Dispose in accordance with governmental environmental regulations. Observe all federal, state, and local environmental regulations.
- **European waste catalogue (EWC):** Not available.

- **Hazardous waste:** Will not occur.

**SECTION 14 – TRANSPORT INFORMATION**

**Note** Please see current shipping paper for most up to date shipping information, including exemptions and special circumstances.

**D.O.T.:** Not regulated.

**IATA.:** Not regulated.

**IMDG:** Not regulated.

**SECTION 15 – OTHER INFORMATION**

**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

**Cholesterol**

CERCLA/SARA Section 313 de minimus %	Not Listed
California Proposition 65	Not Listed
TSCA	Present
EINECS	200-353-2
AICS	Present
Standard for Uniform Scheduling of Medicines and Poisons (SUSMP)	Schedule 4

**1,2-Distearoyl-sn-glycero-3-phosphocholine**

CERCLA/SARA Section 313 de minimus %	Not Listed
California Proposition 65	Not Listed
EINECS	212-440-2

**15.2. Chemical safety assessment**

Chemical Safety Report No information available


- ◆ This product is not classified according to the EU regulations.

**SECTION 16 – OTHER INFORMATION**

- ◆ **Hazardous Material Information System (U.S.A.):**

<b>Health</b>	<b>0</b>
<b>Fire hazard</b>	<b>0</b>
<b>Reactivity</b>	<b>0</b>
<b>Personal protection</b>	

- ◆ **National Fire Protection Association (U.S.A.):**



0 = not significant  
 1 = slight  
 2 = moderate  
 3 = high  
 4 = extreme  
 \* = chronic

- ◆ **Revisions:**

Issue Date: 18-Mar-2026  
 Last Revision Date: 20-Mar-2026  
 Revision note: Name.

- ◆ **Other Comments:**

To the best of our knowledge, the information contain herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

**Restriction of use:** for research use only.